Cavitron Systems, Sterilizable Handpieces and Reservoir Bottles
INFECTION CONTROL PROCEDURES

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The objective of the information provided in this booklet is to supplement published general guidelines for reducing the potential for cross contamination when using a Cavitron® Ultrasonic Scaling and/or Air Polishing Systems, Dispensing System, Cavitron sterilizable handpieces or reservoir bottles during routine dental care. In the event any regulatory agency disagrees with this information, the agency requirements take precedence.

For Infection Control Procedures relating to Cavitron ultrasonic inserts or air-polishing nozzles and inserts, refer to the procedures provided with each instrument.

INFECTION CONTROL PROCEDURES Cavitron Systems

WARNINGS:

- Cavitron system cabinets, power cords, non-sterilizable handpieces, handpiece cables, foot controls, foot control cables, and Sterile Lavage Kit infuser bag are not sterilizable but can be disinfected following the procedures and recommendations listed below.
- Use water based non-immersion type disinfectant solutions. Alcohol based disinfectant solutions may be harmful and may discolor plastic material.

LIMITATIONS ON PROCESSING:

- Do not spray disinfectant solutions directly on system surfaces.
- Do not use disinfectant solution on sterilizable handpieces. Refer to sterilization procedures for sterilizable handpieces.

INSTRUCTIONS:***Point of Use:***

Remove excess soil with disposable cloth or paper. Discard cloth after use.

Containment and Transportation:

Handle with care.

Cleaning: manual

Generously spray disinfectant solution on a clean cloth. Wipe the outer surfaces of the cabinet, power cords, non-sterilizable handpieces, handpiece cables, air and water lines, foot control, foot control cables, and infuser bag. Discard used cloth. Wipe dry with a clean cloth.

Disinfection:

Generously spray disinfectant solution on a clean cloth. Wipe the outer surfaces of the cabinet, power cords, non-sterilizable handpieces, handpiece cables, air and water lines, foot control, foot control cables, and infuser bag. Discard used cloth. Allow disinfectant to air dry.

Drying:

When cleaning, wipe surfaces dry with a clean cloth. To achieve disinfection, allow surfaces to air dry.

Maintenance, inspection and testing:

- Visually inspect to ensure all contamination is removed.
- Visually inspect all cables and power cord for tears.

- Visually inspect handpiece cable o-rings for wear. Replace worn o-rings (refer to applicable Directions for Use of system). Lubricate o-rings with water before use.

Storage:

Do not place on or next to a radiator or other heat source. System is portable but must be handled with care.

INFECTION CONTROL PROCEDURES

Sterilizable Handpieces, Handpiece Adapters, and Reservoir Bottles

WARNINGS:

- Instruments are not sterile upon receipt and must be sterilized prior to use in accordance with the following instructions.
- The colorant used on handpieces with soft grips may bleed during autoclaving. This will not effect performance.
- Do not disinfect or autoclave the Steri-Mate Handpiece or Steri-Mate Handpiece Adapter while assembled.

LIMITATIONS ON PROCESSING:

- Repeated cleaning and sterilization cycles have minimum effect on these instruments. End of life in normally determined by wear and damage due to use.
- Reservoir bottle storage caps are not steam autoclave sterilizable.
- Do not use chemical disinfectants prior to sterilization or rapid deterioration of the material may occur.

- Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.
- Do not immerse handpieces or handpiece adapter in ultrasonic bath.

INSTRUCTIONS:

Point of Use:

Remove excess soil with disposable cloth or paper.

Containment and Transportation:

Protect instruments from damage. Instruments should be reprocessed as soon as is reasonably practical following use. Instruments left standing wet may stain or corrode.

Preparation for decontamination:

Remove excess powder from the air polishing or air polishing/scaling handpieces by inserting the applicable cleaning wire(s) provided with your system or handpiece into the powder ports and moving the wire back and forth several times. Remove wire for cleaning and sterilization.

Cleaning: automated

Reservoir bottles can be washed in a dishwasher.

Cleaning: manual

- Rinse all instruments under running water to remove any gross debris.
- Wipe sterilizable handpieces, reservoir bottles and storage caps, blue interface caps and sip tubes

with a cloth saturated with multi-purpose cleaner and rinse with running water. Discard used cloth. Dry with a clean cloth.

Disinfection:

To disinfect the storage caps, generously spray disinfectant solution on a clean cloth and wipe the caps thoroughly. Discard used cloth. Allow disinfectant to air dry. Storage caps are not steam autoclave sterilizable.

Disinfection of the sterilizable handpieces, handpiece adapters, reservoir bottles, blue interface caps and sip tubes is not necessary prior to steam autoclave sterilization.

Packaging:

Place instrument in a paper or paper/plastic steam sterilization pouch. Instruments may be loaded into dedicated instrument trays.

Sterilization:***(Handpieces, handpiece adapters, Reservoir Bottles, and blue interface caps with sip tubes)***

Use a steam autoclave. Place bagged instruments into the steam autoclave, paper side up when using a paper/plastic pouch. Disassemble reservoir bottles and interface caps before placing in pouch. After warm-up is completed, operate at a sterilizing temperature of 273°F/134°C to 279°F/137°C and pressure of 31psi/216kPa for 3 to 12 minutes, except when sterilizing the GRAY Steri-Mate handpiece, the cycle time must be 4-12 minutes.

Alternate method: Place non-bagged instruments into the steam autoclave. After warm-up is completed, operate at a sterilizing temperature and pressure of 273°F/134°C to 279°F/137°C and pressure of 31psi/216kPa for 3 to 4 minutes, except when sterilizing the GRAY Steri-Mate handpiece, the cycle time must be 4-5 minutes. Instruments should be used immediately.

Drying:

To dry, use the drying cycle of the autoclave. Set cycle for 20-30 minutes.

Maintenance, inspection and testing:

- Visually inspect to ensure all contamination is removed.
- Check for distortion, damage and wear. Discard damaged, worn or corroded instruments.
- Inspect o-rings on the handpiece, handpiece adapter and handpiece cable for damage or wear. Replace if necessary.
- Ensure all o-rings are in place.

Storage:

To maintain sterility, instruments should remain bagged until ready for use.

MANUFACTURER CONTACT:

In the United States, contact DENTSPLY Professional Customer Service Technical Support at 800-989-8826. For areas outside the United States, contact your local Dentsply Professional Representative.

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DENTSPLY DeTrey GmbH