ChemFil® Superior

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ChemFil® Superior

Fast-Setting Glass Ionomer Restorative Material

ChemFil® Tooth Cleanser

Cleanser for Cavity Conditioning

ChemVarnish™

Protective Varnish for Glass Ionomer Restorations

CAUTION: For dental use only.
USA: Rx only.

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	5
4 Hygiene	7
5 Lot number and expiration date	7

1 Product description

ChemFil® Superior restorative is a fast-setting glass ionomer restorative material consisting of a blend of alumino silicate glass and polyacrylic acid. The powder is mixed with distilled water to produce a filling material which adheres to dentin and enamel producing tightly sealed, strong and esthetic restorations.

ChemFil® Superior restorative is available in seven shades:

- L – Light, similar to VITA®¹ shade A1
- LY – Light-Yellow, similar to VITA® shade C1
- LYG – Light-Yellow-Grey, similar to VITA® shade A2
- LG – Light-Grey, similar to VITA® shade B2
- DG – Dark-Grey, similar to VITA® shade D4
- DY – Dark-Yellow, similar to VITA® shade A3
- GB – Grey-Brown, similar to VITA® shade A4

ChemFil® Tooth Cleanser may be used for conditioning of the cavity surface to enhance cavity cleaning and therefore adhesion.

ChemVarnish™ is recommended to protect glass ionomer restorations from moisture contamination and dehydration during or immediately after setting.

¹ Not a registered trademark of Dentsply Sirona, Inc.

1.1 Delivery forms

ChemFil® Superior:

- Powder/liquid for hand mixing.

ChemFil® Tooth Cleanser:

- 6 ml bottle.

ChemVarnish™:

- 12.5 ml bottle.

1.2 Composition

ChemFil® Superior:

- Calcium-alumino-sodium-fluoro-phosphor-silicate glass
- Polycarboxylic acid
- Tartaric acid
- Iron oxide pigments

ChemFil® Tooth Cleanser:

- Polycarboxylic acid
- Water

ChemVarnish™:

- Polyethyl methacrylate
- Ethyl acetate

1.3 Indications

Suitable for:

- Restoration of Class I and II cavities of deciduous teeth.
- Restoration of Class III cavities.
- Restoration of Class V lesions and cavities.
- Temporary and semi-permanent restorations.
- Fissure fillings (minimal Class I cavities).

1.4 Contraindications

- Direct or indirect application to dental pulp (direct or indirect pulp-capping).
- Permanent restoration of occlusal stress-bearing areas.

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

ChemVarnish™ contains ethyl acetate and may be irritating to skin, eyes and oral mucosa.

ChemFil® Superior and ChemFil® Tooth Cleanser are acidic and may cause burns.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse immediately with plenty of water and seek medical attention.

- **Avoid skin contact** to prevent irritation. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and wash thoroughly with soap and water.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- To prevent bottles, powder scoop and mixing pad from exposure to spatter or spray of body fluids or contaminated hands it is recommended that the bottles, powder scoop and mixing pad are handled offside the dental unit with clean/disinfected gloves.
- The shade guide cannot be reprocessed. To prevent the shade guide from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the shade guide is handled with clean/disinfected gloves.
- Applicator Tips are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- The desiccant crystals integrated in the cap and/or drying capsule protect the powder from absorbing moisture. The color of the desiccant crystals should be orange. The desiccant crystals are no longer effective after the color has turned from orange to colorless.
- Tightly close bottle immediately after use.
- Increased temperature and/or humidity reduce available working time and accelerate setting time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.
- The product should be mixed in proportions of two level scoops powder to two drops water for optimum performance. The mixed material should be homogeneous and streak free prior to application. Variations may affect the strength and durability of the material.
- Use in a well-ventilated area. Avoid ingestion or prolonged inhalation.
- ChemVarnish™ is flammable: Keep away from sources of ignition.
- Do not empty into drains.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store with original caps tightly closed in a well-ventilated place at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 Shade selection

ChemFil® Superior restorative offers seven shade options:

- L – Light, similar to VITA® shade A1
- LY – Light-Yellow, similar to VITA® shade C1
- LYG – Light-Yellow-Grey, similar to VITA® shade A2
- LG – Light-Grey, similar to VITA® shade B2
- DG – Dark-Grey, similar to VITA® shade D4
- DY – Dark-Yellow, similar to VITA® shade A3
- GB – Grey-Brown, similar to VITA® shade A4

Shade selection is made with the shade guide supplied. As the translucency and color of glass ionomer cements change during the first hours after placement, the shade guide represents the final shade of the restoration.

3.2 Cavity preparation

1. In all cavity classes cavity preparation may be kept to the minimum required for caries removal.
2. Preparation of mechanical retention is not needed except when restoration is liable to be subjected to high displacement forces.
3. Rinse surface with water spray and carefully remove excess rinsing water with air spray. Do not desiccate the tooth substance – cavity surface should be moist.
4. Use cotton rolls or rubber dam to isolate the cavity from contamination.

3.3 Cleaning

1. Clean uninstrumented enamel and dentin with a rubber cup and pumice or a cleaning paste such as Nupro® prophylaxis paste.
2. Wash thoroughly with water spray and air dry.
3. Clean freshly instrumented enamel and dentin with water spray, then air dry.

3.4 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see complete Instructions for Use) leaving the rest of the cavity surface free for bonding.

3.5 Conditioning tooth surface

For surface conditioning, ChemFil® Tooth Cleanser may be used.

1. Place an appropriate amount of ChemFil® Tooth Cleanser onto the pad provided or on a glass slab.
2. Apply ChemFil® Tooth Cleanser to the cavity surface using an Applicator Tip.
3. Allow ChemFil® Tooth Cleanser to remain for 15 seconds.
4. Rinse surface with water spray and carefully remove excess rinsing water with air spray. Do not desiccate the tooth substance – cavity surface should be moist.

3.6 Placement of matrix

1. In Class II cavities place a matrix (e.g. AutoMatrix® matrix system or Palodent® Plus matrix system) and wedge. Burnishing of the matrix band will improve contact and contour.
2. Contact surface of steel matrix can be coated with a water soluble or petroleum lubricating jelly.
3. In Class V cavities cervical matrices may be used.

3.7 Dosage and mixing

The mixing ratio is 2 scoops powder : 2 drops liquid².

1. Invert powder bottle to fluff the powder.
2. Overfill the scoop with powder and strike off the excess with the bottle insert leaving a flat surface level with the edge of the scoop.

NOTICE

Premature setting of material.

Reduced working time.

1. Keep powder and liquid separated until ready for mixing.
3. Place the powder onto the pad provided or on a glass slab³. Special spatulas are not required.
4. Divide the powder into two halves.
5. Fill water dispenser with distilled or de-ionized water.
6. Dispense two drops of water by holding the bottle vertically and squeezing gently.
7. Mix the first part with the water in five seconds or less, then add the second part and mix for about ten seconds – **maximum mixing time 20 seconds**. The consistency of the mixed material should resemble that of a composite resin.

3.8 Working and setting time⁴

- From start of mixing, the remaining working time is approximately 2 minutes.
- After mixing, the setting time is 2 to 3 minutes.

3.9 Application

1. Avoid contamination of the filling with saliva or water during insertion.
2. Pack the mixed cement firmly into place and contour with a suitable matrix.
3. Keep the matrix in place for at least 4 minutes.
4. Immediately upon removal of the matrix, coat the restoration with ChemVarnish™.
5. Place an appropriate amount of ChemVarnish™ onto the pad provided or on a glass slab.
6. Apply ChemVarnish™ to the restoration surface using an Applicator Tip.
7. Dry varnish layer using the air syringe.
8. To ensure a continuous film apply a second layer of varnish and dry again.

3.10 Finishing

Between 7 minutes and 24 hours after placement:

Finishing and polishing is best carried out at the next appointment. However, the restoration is sufficiently hard for water-free finishing approximately 7 minutes after placement.

1. Finishing and polishing instruments should be lubricated with Vaseline. Do not use water spray.
2. Before dismissing the patient the restoration should be protected from dehydration and moisture by painting the surface with ChemVarnish™.

If finished more than 24 hours after placement:

1. Water spray may be used during finishing.
2. Restoration need not be protected with ChemVarnish™.
3. Remove any large excess of material with finishing burs and complete finishing with finishing discs or strips.
4. Finally buff the restoration to a high lustre with polishing cups and paste. The Enhance® Finishing and Polishing System is recommended.

² The same powder : liquid ratio (7.4 g : 1 g) is proposed for testing purposes according to ISO 9917.

³ Mixing and application of the material should be done at normal room temperature. Higher temperatures shorten the working and setting times of the cement. Under warm conditions or if a long working time is required, it is therefore recommended that the glass slab be cooled. Testing to ISO 9917 is carried out at 23 °C ± 1 °C and a relative humidity of 50 ± 10%.

⁴ Net setting time determined according to ISO 9917.

3.11 Cleaning instruments

Glass ionomers adhere to instruments and glass slabs.

1. Immediately after use, either clean instruments with water or drop in disinfectant or water to facilitate subsequent cleaning.

4 Hygiene

Cross-contamination.

Infection.

1. The shade guide cannot be reprocessed.
2. To prevent the shade guide from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the shade guide is handled with clean/ disinfected gloves. Do not reuse shade guide if contaminated.
3. Dispose contaminated shade guide in accordance with local regulations.

Cross-contamination.

Infection.

1. Do not reuse single use products. Dispose in accordance with local regulations.

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2014-11-28

[These Instructions for Use are based on Master Version 04]

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Dentsply DeTrey GmbH