# Dyract<sup>®</sup> flow

## Flowable Compomer Restorative

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CAUTION: For dental use only.

USA: Rx only

Content	Page
1 <b>Product description</b> .....	2
2 <b>Safety notes</b> .....	3
3 <b>Step-by-step instructions</b> .....	5
4 <b>Hygiene</b> .....	7
5 <b>Lot number and expiration date</b> .....	7

### 1 Product description

**Dyract<sup>®</sup> flow** restorative is a compomer restorative material with flow characteristics which make it ideal for small cavities in anterior and posterior teeth and for cavity lining. **Dyract<sup>®</sup> flow** restorative adapts to the cavity walls without the use of hand instruments.

**Dyract<sup>®</sup> flow** restorative is available in 5 VITA<sup>®</sup><sup>1</sup> shades in a regular opacity (A2, A3, A4, B1) and an opaque (O-A3) shade.

#### 1.1 Delivery forms

- Small syringes containing 1.8 g (1 ml) restorative to be used in combination with disposable steel canulas (applicator needles, dispensing tips) for direct intraoral application.

#### 1.2 Composition

- Strontium-alumino-fluoro-phosphor-silicate glass
- Highly dispersed silicon dioxide
- Ammonium salt of phosphoric acid modified methacrylate resin
- Carboxylic acid modified methacrylate resin
- Diethylene glycol dimethacrylate (DGDMA)
- Camphorquinone
- Ethyl-4(dimethylamino)benzoate
- 2-Hydroxymethoxybenzophenone
- Butylated hydroxy toluene (BHT)
- Ironoxide pigments
- Titanium dioxide

<sup>1</sup> Not a registered trademark of Dentsply Sirona, Inc.

### 1.3 Indications

- Restoration of minimally invasive cavity preparations, pits and fissures, including air abrasion and tunnel preparations.
- Small restorations without antagonist contact in anterior and posterior teeth.
- Restoration of minimal, shallow Class V preparations, including incipient Class V abfraction and erosion lesions.
- Base/liner under composite, high-density composite and compomer Class I & Class II restorations.
- Undercut blockout and filling of small defects in indirect restoration preparations.
- Small margin defect repair in otherwise sound restorations.

### 1.4 Contraindications

- Use with patients who have a history of severe allergic reaction to methacrylate resins.

### 1.5 Compatible adhesives

The material is chemically compatible with methacrylate-based dentin/enamel adhesives including Dentsply Sirona adhesives designed for use with visible light cured composite restoratives (see complete Instructions for Use of selected adhesive).

Xeno<sup>®</sup>, self etching adhesives and Prime&Bond<sup>®</sup> NT, etch and rinse adhesive are especially recommended.

## 2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

#### Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

### 2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

## 2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Only use the canulas supplied by Dentsply Sirona for placement of the material.
- Canulas are used for direct intraoral application of restorative. Canulas are for single use only, discard in order to prevent cross-contamination. Canula may clog as restorative may set inside.
- Tightly close syringes immediately after use. Do not allow canula to remain as a syringe cap. Replace original cap tightly after each use.
- The material should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may cause material and not fully seated canulas to eject from the syringe.
- Interactions:
  - Do not use eugenol- and hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening and cause softening of the product.

## 2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

## 2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in sealed packs at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Keep out of direct sunlight and protect from moisture.
- Do not freeze.
- Do not use after expiration date.

Humidity can adversely affect the properties of unsealed syringes. For unsealed syringes the following is recommended:

- Store in a dry environment (< 70% relative humidity).
- Use unsealed syringes within 5 months.

### 3 Step-by-step instructions

#### 3.1 Shade selection

Shade selection is made whilst the teeth are hydrated and prior to the restorative procedure.

1. Remove any extraneous plaque or surface stain.
2. Use the VITAPAN®<sup>2</sup> classical shade guide. Use the central part of the respective VITA® tooth for shade selection.

#### 3.2 Cavity preparation

1. Clean the tooth surface with a fluoride free prophylaxis paste (e.g. Nupro® prophylaxis paste).
2. Prepare cavity (unless not necessary, e.g. cervical lesion).
3. Use adequate isolation such as rubber dam.
4. Rinse surface with water spray and carefully remove rinsing water. Do not desiccate the tooth structure.

#### 3.3 Matrix placement and wedging (Class II restorations)

Place a matrix (e.g. AutoMatrix® matrix system or Palodent® V3 Sectional matrix system) and wedge. Burnishing of the matrix band will improve contact and contour.

#### 3.4 Pulp protection, tooth conditioning/dentin pre-treatment, adhesive application

Refer to adhesive manufacturer's Instructions for Use for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of the material.

#### 3.5 Application

Apply material in increments immediately after the application of the adhesive. The material is designed to be cured in increments up to a 2 mm depth/thickness.

1. Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto a mixing pad.
2. Attach disposable steel canula<sup>3</sup> to the end of the syringe. Turn canula clockwise until fully seated. Tug on canula to assure locking into syringe collar.
3. Dispensing tip may be bent for easy access. Gentle curving of the tip is preferable to an acute bend. Bending may occlude the canula. Be sure to check that this has not occurred by extruding a small amount through the tip.
4. Material should flow freely with GENTLE PRESSURE. DO NOT USE EXCESSIVE FORCE. If not, remove syringe from patient field and check for obstruction.
5. Apply material through the canula to surfaces to be restored.

#### **Injury due to excessive force.**

1. Apply slow and steady pressure on the syringe.
2. Do not use excessive force. Excessive pressure may cause material and not fully seated canulas to eject from the syringe.

<sup>2</sup> Not a registered trademark of Dentsply Sirona, Inc.

<sup>3</sup> Use the disposable steel canulas included in the Dyract® flow package, ReOrder No.: 606.04.490.

### 3.5.1 Instructions for application as cavity liner

For most cavity lining applications, a thickness of 1 mm or less is recommended. With the exception of Class II proximal box (gingival) margins, cavity margins/conditioned enamel are usually closed with selected restorative material. Lining is generally confined to internal conditioned dentin surfaces. Avoid lifting the tip while dispensing to minimize air entrapment. Any visible air bubbles should be pierced with a clean, sharp explorer prior to curing. Technique Tip: Use of a clean instrument such as a sharp explorer may be used to evenly spread a thin layer of dispensed Dyract® flow restorative material over prepared cavity (dentin) surfaces prior to curing.

Proceed immediately to complete the restoration per restorative material manufacturer's Instructions for Use. For adhesive restorations, do not disturb or allow contamination of cured liner surface. If contamination occurs, follow adhesive manufacturer's Instructions for Use for re-application of adhesive.

Technique Tip: If surfaces remain undisturbed, additional application of adhesive between increments is neither necessary nor recommended.

### 3.6 Light curing

1. Light cure each area of the restoration surface with a suitable curing light designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 500 mW/cm<sup>2</sup> exposure for at least 20 seconds. Refer to curing light manufacturer's Instructions for Use for compatibility and curing recommendations.
2. Cure each layer according to the table below.

#### Inadequate polymerization due to insufficient curing.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output.

#### Curing Light Output $\geq$ 800 mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
A2, B1	2 mm	10 sec
A3, A4	2 mm	20 sec
O-A3	2 mm	30 sec

#### Curing Light Output 500-800 mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
A2, A3, B1	2 mm	20 sec
A4, O-A3	2 mm	40 sec

### 3.7 Finishing and polishing

1. Contour the restoration using finishing burs or diamonds.
2. Use Enhance® finishing and polishing instruments and interproximal finishing strips.

For contouring, finishing, and/or polishing, follow the manufacturer's Instructions for Use.

## 4 Hygiene

### 4.1 Protecting and disinfecting the syringes

1. Tightly close syringes immediately after use. Do not allow canula to remain as a syringe cap. Replace original cap tightly after each use.
2. To prevent syringes from exposure to spatter or spray of bodily fluids, contaminated hands or oral tissues, use of a protective barrier is recommended.
3. Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.

#### **Damage to the label by repeated disinfection.**

Avoid contamination of syringes.

### 4.2 Steel canulas – single use

#### **Cross-contamination of steel canulas.**

Do not reuse steel canulas.

Properly dispose steel canulas in accordance with local regulations.

## 5 Lot number ( ) and expiration date ( )

1. Do not use after expiration date.  
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
  - Reorder number
  - Lot number
  - Expiration date

© Dentsply Sirona 2010-02-25

[These Instructions for Use are based on Master Version 01]

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Dentsply DeTrey GmbH