Enhance[®] Multi

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Enhance[®] Multi

Multiple Use Finishing & Polishing System

CAUTION: For dental use only.

USA: Rx only

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	5
4 Hygiene	5
5 Lot number and expiration date	7

1 Product description

The **Enhance[®] Multi** Finishers and Polishers are multiple use, diamond impregnated, silicone Discs, Cups and Points designed for preparing and polishing the surface of composite and compomer restorations in a two step technique.

1.1 Indications

- The dark grey Enhance[®] Multi Finishers are indicated for contouring and pre-polishing of composite and compomer restorations.
- The light grey Enhance[®] Multi Polishers are indicated for polishing of composite and compomer restorations.

1.2 Contraindications

- Not known.

1.3 Delivery forms

- Assortment of Enhance[®] Multi Finishers, Polishers and accessories plus organizer tray.
- Refill packages containing 5 instruments of the same type each.

1.4 Composition

Active part:

- Diamond grit
- Synthetic rubber
- Pigments

Shaft:

- Stainless steel

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use. Improper use of these instruments leads to bad quality results and more risks. This is why they must be used only by qualified persons.

2.1 Warnings

- Enhance® Multi Finishers and Polishers are designed for use with water. As with any rotary instrument, heat will build up with prolonged contact. Use at a speed between 5,000 and 10,000 rpm maximum with light intermittent pressure. Never exceed the maximum rotation speed. Higher speed may result in heat buildup and may endanger the stability of the instruments causing danger to user and patient. Heat buildup may cause pulp inflammation.
- In case of non-observance of the maximum rotation speed allowed, the polisher will tend to vibrate. Such vibrations can deform the polisher or its shaft and cause them to break. In the latter case, the user, the patient and other persons can be injured.
- Care should be taken to have good control on any rotary instrument to protect patients from injury. Contact with soft tissues (skin, gingiva, mucosa) may result in injuries. If injury occurs, ensure appropriate medical care.
- Dental professionals should wear breathing protection (surgical mask). We recommend a dust aspiration system.
- Dental professionals and patients should wear protective eyewear. In case of improper use or material failure, the mandrel, shaft or workpiece can break and could become a dangerous flying object.
- Excessive pressure on the polisher can destroy it.
- Excessive pressure results in more heat being produced.
- Excessive pressure can cause early wear of the polisher.

2.2 Precautions

This product is intended to be used as specifically outlined in the Instructions for Use.

Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.

- Rotating instruments have to be cleaned, disinfected, dried and sterilized prior to the first application on the patient and immediately after every use.
- These instruments must only be used with turbines, handpieces and contra-angles that are technically and hygienically perfect, meaning that they should be well maintained and correctly cleaned. The turbines and contra-angles used for these instruments must ensure precise and regular rotation.
- Insert the instruments as deeply as possible without excessive use of force into the chuck and check to ensure that they are secure. Before applying the instrument on the surface to be treated, bring it to the rotation speed needed.
- In order to avoid deformations, use circular movements with the instrument. Do not bend the instrument or use it as a lever because this would lead to increased risk of breaking the instrument.
- Immediately remove any instruments that are damaged, bent or no longer run concentrically.
- Unmounted polishers must be centred after shaftmounting in order to avoid vibrations during use and only high quality mandrels should be applied. Using mandrels of inferior quality could break them and cause injuries.
- Observing rotation speeds ensures optimum results. Non-observance of the maximum rotation speed allowed increases safety risks.
- Always avoid excessive pressure on the instrument because this will cause overheating, which could damage the pulp. Additionally, in case of very excessive pressure, the instrument can break and, as a result, persons can be injured.
- In order to avoid the overheating of a tooth, sufficient water cooling must be ensured (50 ml/min). Insufficient water cooling can lead to irreversible damage to the tooth and its surrounding tissues.
- Moistening of the surfaces to be prepared must not be impaired by an incorrect suction technique: Deflection of the coolant jet can cause heat related damage.
- Protect adjacent restorations from accidental contact with rotating disc, cup or point.

- Excessive speeds while using instruments with overly large heads endanger the stability of the instruments.
- To maintain the longevity of these instruments, careful preparation with low contact pressure is important: the recommended contact pressure range is 0.3-2 N (30-200 p).
- Avoid tilting, levering, or blocking the instruments up by using too much contact pressure (increased danger of fracture).
- The use of inappropriate disinfectants and/or prolonged disinfectant times can result in corrosion (increased danger of breakage) and discoloration of the instruments. Do not use alkaline or aldehyde containing disinfectants. Do not place in alcohol.
- Use in a well-ventilated area. Avoid inhaling airborne particulates.

2.3 Adverse reactions

- Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury. (See 2.1 Warnings and 2.2 Precautions).
- Product misuse may cause excessive heat. (See 2.1 Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Dry storage is recommended.
- Store at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Store the Enhance® Multi Finishers and Polishers protected against contaminants and against chemicals, acids, heat, mechanical stress and extreme temperature variations.
- Do not use after expiry date.

3 Step-by-step instructions

Enhance® Multi Finishers and Polishers are only to be used with a RA shank handpiece.

1. Place restoration or composite veneer as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using finishing burs.
2. Insert an Enhance® Multi Finisher (dark grey) into a conventional speed contraangled hand piece and continue finishing. Use with water. The aggressiveness of the Enhance® Multi Finishers is controlled by the pressure applied to the surface of the composite. Use moderate to light, intermittent pressure in a buffing motion.
Frequently raise instrument from surface contact.
3. Once final contours are achieved, proceed to polish. Insert an Enhance® Multi Polisher (light grey) and continue polish for final high surface luster.

Note: Enhance® Multi Finishing and Polishing Points are ideal for posterior occlusal and concave lingual surfaces. Enhance® Multi Finishing and Polishing Cups are well suited to cervical and proximal line angle areas.

4 Hygiene

Disposable products delivered unsterile, marked with the symbol (⊗), may only run through the validated sterilization cycle ONCE before initial use.

Concerns all rotating polishing and grinding instruments that are classified as semi-critical according to the RKI-guideline.

- All instruments are delivered unsterile and must be run through the indicated cycle before and after each use.
- The label on the cleaning and/or disinfecting solution must specifically state "suitable for rubber polishers or synthetics/silicones", because if the label only states "rotary instruments", the solution may not be suitable for polishers.

4.1 Warning notices

- Strong acids and strong bases may oxidize the stainless steel shaft.
- Rinse the polisher with distilled water after the treatment with cleaning and disinfectant solutions.
- Avoid temperatures >150 °C.
- Ultrasonic bath must not exceed temperatures of 42 °C because of the possible coagulation of albumen.

4.2 Reprocessing instructions

Cleaning preparation	Pre-clean under running water with a brush (plastic) directly after use.
Manual cleaning	Under running water with a brush (plastic bristles). Machine assisted: 1. Ultrasonic assisted cleaning with a suitable cleaning agent and disinfectant. 2. Rinse the instruments afterwards under running water.
or	
Automatic cleaning	Thermal disinfectant: Manufacturer's specification according to DIN EN ISO 15883. Cleaning program as indicated by the manufacturer in the operating instructions.
Disinfection	Use a solution classified as suitable for rubber and silicone polishers and synthetics by the disinfectant manufacturer. Exposure times and concentrations recommended by the manufacturer should be adhered to.

Rinsing	Rinse the polishers with distilled water after treatment with cleaning and disinfectant solutions.
Drying	Dry with fresh, clean, lint-free cellulose tissues.
Maintenance	Visual check of all instruments with optical magnification (5-10 fold).
Inspection	No residues > continue to sterilization . Visible residues > repeat cleaning. Reject and dispose of instruments in the event of discernible defects.
Sterilization	For all instruments that need to be sterilized in accordance with EN ISO 17664 and all national valid legal requirements. Steam sterilization: Appliance according to EN 13060, validated procedure. Category: S- or B-sterilizer. <ul style="list-style-type: none"> • Holding time: Full cycle 5 min • Sterilization temperature: 134 °C • Drying time: 10 min • Threshold values of contents for feed-water and steam condensates • Loading of sterilizer according to manufacturer's instructions • Follow manufacturer's operating instructions
Storage	Store instruments packed and protected from recontamination in proven suitable sterile packaging, cassettes or retainers.
Additional information	Repetitive reprocessing can change both the look and feel of the instruments slightly, but does not interfere with the instrument's function.

4.3 Validation conditions

Manual cleaning:

- Cleaning and disinfectant solution Dürr Dental ID 212
- Concentration: 2%
- Holding time: 5 min

Ultrasonic assisted cleaning:

- Dürr Dental ID 212/Readymade solution; Holding time according to manufacturer
- Type of appliance: Ultrasonic
- Concentration: 2%
- Holding time: 2 min

Automatic cleaning:

- Cleaning and disinfecting machine Miele G 7883
- Program: SPECIAL 93 °C – 10'
- Holding time: 10 min
- Detergent: Dr. Weigert – neodisher MediClean Dental
- Rinsing: 3 min at 75 °C with neodisher Z Dental

Sterilization:

- Holding time full cycle: 5 min
- Sterilization temperature: 134 °C
- Drying time: 10 min
- Threshold values of contents for feed-water and steam condensates

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2017-05-11

[These Instructions for Use are based on Master Version 04]

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Dentsply DeTrey GmbH